fr	NOTICE D'UTILISATION – PROTECTEUR OCULAIRE	2
en	INSTRUCTION FOR USE – OCULAR SHIELD	2
de	GEBRAUCHSANWEISUNG – OKULARSCHUTZKAPPE	3
it	ISTRUZIONI PER L'USO – SCUDO OCULARE	4
es	INSTRUCCIONES DE USO – PROTECTOR OCULAR	5
pt	INSTRUÇÕES DE USO – PROTETOR OCULAR	6
nl	GEBRUIKSAANWIJZING – OOGBESCHERMER	7
da	BRUGSVEJLEDNING – ØJENSKJOLD	8
fi	KÄYTTÖOHJEET – SILMÄSUOJUS	9
sv	BRUKSANVISNING – OKULÄR SKÖLD	10
el	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ – ΕΝΔΟΦΘΑΛΜΙΟ ΚΑΛΥΜΜΑ	11
et	KASUTUSJUHEND – SILMAKILP	12
lv	LIETOSANAS INSTRUKCIJA – ACS AIZSARGS	13
lt	NAUDOJIMO INSTRUKCIJOS – AKIES APSAUGOS PRIEMONĖ	14
pl	INSTRUKCJA UŻYTKOWANIA – OSŁONA OCZU	15
sl	NAVODILA ZA UPORABO – ŠČITNIK ZA OČI	16
sk	NAVOD NA POUZIVANIE – OČNÝ ŠTÍT	17
cs	NAVOD K POUZITI – OČNÍ ŠTÍT	18
ro	INSTRUCȚIUNI DE UTILIZARE – PROTECTOR OCULAR	19
bg	ИНСТРУКЦИИ ЗА УПОТРЕБА – ОЧЕН ПРОТЕКТОР	20
hr	UPUTE ZA UPOTREBU – OČNI ŠTIT	21
hu	HASZNALATI UTASITAS – SZEMVÉDŐ	22
no	BRUKSANVISNING – ØYESKJOLD	23
tr	KULLANIM – OKÜLER KALKAN	24
zh-tw	使用说明书 – 眼罩	25
Illustrations.....		26
Symbols.....		27
IFU-2138 311 18 D		

FCI S.A.S.
20-22, rue Louis Armand
75015 PARIS – FRANCE

: + 33 (0) 1 53 98 98 98

: + 33 (0) 1 53 98 98 99

@ : fci@fci.fr
 : www.fci.fr

Distributed in U.S.A. by:

FCI Ophthalmics Inc.
30 Corporate Park Drive
Suite # 310/320
Pembroke, MA 02359

: 800-932-4202

: 781-826-9062

@ : info@fci-ophthalmics.com
 : www.fci-ophthalmics.com

Rx Only

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Product scope of this IFU:

- S6.2301: OCULAR SHIELD

OCULAR SHIELD is a protective egg-shaped shield (maximal diameter: 27mm) made of PMMA (Polymethylmetacrylate), with a prehension pin which is used to place and remove the shield easily.

The device is sterilized using Gamma radiation in its final simple blister packaging. A red indicator shows that the product followed a sterilization cycle validated by FCI. The device is supplied by box of 1. The device is for single use only and is not resterilizable.

Intended Purpose

OCULAR SHIELD is designed to protect the eyeball during conventional eyelid surgery or lacrimal duct surgery (DCR).

Indications

Medical conditions: patients undergoing a periorbital surgery

Intended patient population: adults

Body parts in contact with the device: internal surface of the eyelid, cornea and sclera

Intended users: ophthalmic surgeons

Intended use environment: operating room

Contraindications

The OCULAR SHIELD should not be used during a surgery involving a laser.

Warnings

The device is for single use and must not be resterilized. Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety.

Check packaging integrity before use, and do not use the device in the event of packaging being damaged, unintentionally opened before use, or if it was exposed to storage or transport conditions outside of those specified.

Do not use the product if the indicator is not red, as it could endanger the patient. Do not use after the expiration date shown on the package.

The device is intended to be removed at the end of the intervention.

In case the user identifies faulty labelling information, it shall be relayed to FCI S.A.S. and the device should not be used if the patient's safety is threatened.

Precautions for use

the device must be stored at temperatures between 5°C/41°F and 35°C/95°F and humidity between 30% and 75%, away from sunlight.

The device may be exposed to temperature variations between -15°C/5°F and 50°C/122°F and humidity variations between 15% and 90% during transport phase for 10 days.

Remove the device from its packaging under aseptic conditions.

Adverse effects

As in any type of surgery, there are risks linked to the material or to developments of the initial pathology. Potential complications associated with the use of the device include, but are not limited to the following:

- Corneal infection
- Sub conjunctival hemorrhage
- Corneal abrasion
- Red eyes and photophobia
- Discomfort

Any complication that occurs in relation with the device should be notified to FCI S.A.S. and the competent authority of the Member State in which the user and/or the patient is established, or, for users and/or patients not subject to EU requirements, to FCI S.A.S. and to the applicable regulatory authority in conformity with the requirements of that regulatory authority

Instructions for use:

Use the OCULAR SHIELD prehension pin to place and remove it easily.

Performances, mode of action and clinical benefits:

OCULAR SHIELD forms a mechanical barrier between the eyelids and the cornea, which allows the surgeon to protect the eyeball of the patient during the surgery.

Compatibility with other devices

FCI S.A.S. does not claim the compatibility of OCULAR SHIELD with other devices.

Information to be communicated to the patient

Any serious incident that occurs in relation with the device should be notified to FCI S.A.S. and the competent authority of the Member State in which the patient is established.

The user must inform the patient of the contraindications and potential complications related to the device that are set out in the present document.

Waste disposal

With household waste or recycle	With hospital waste (contaminated products)
Outer box, instructions for use	Ocular shield
Blister	

The disposal of this device does not present any physical risk

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